



## DOSAGE GUIDE

### NuWell® 410 CHLORINE ENHANCER

| Nominal Well Size | Well Volume in<br>(Gal/Ft) | Amount of NuWell<br>410 (Quarts/Ft) | Amount of Chlorine Product |                        |                       |
|-------------------|----------------------------|-------------------------------------|----------------------------|------------------------|-----------------------|
|                   |                            |                                     | CaHypocl-65%<br>(Lbs/Ft)   | SodHyp-12%<br>(Gal/Ft) | SodHyp-5%<br>(Gal/Ft) |
| 2                 | 0.16                       | 0.0007                              | 0.0004                     | 0.0003                 | 0.0006                |
| 3                 | 0.37                       | 0.0015                              | 0.0010                     | 0.0006                 | 0.0014                |
| 4                 | 0.65                       | 0.0026                              | 0.0017                     | 0.0011                 | 0.0025                |
| 5                 | 1.02                       | 0.0041                              | 0.0027                     | 0.0017                 | 0.0039                |
| 6                 | 1.47                       | 0.0059                              | 0.0038                     | 0.0025                 | 0.0056                |
| 8                 | 2.62                       | 0.010                               | 0.007                      | 0.004                  | 0.010                 |
| 10                | 4.09                       | 0.016                               | 0.011                      | 0.007                  | 0.016                 |
| 12                | 5.89                       | 0.024                               | 0.015                      | 0.010                  | 0.022                 |
| 14                | 8.02                       | 0.032                               | 0.021                      | 0.014                  | 0.030                 |
| 16                | 10.47                      | 0.042                               | 0.027                      | 0.018                  | 0.040                 |
| 18                | 13.25                      | 0.053                               | 0.034                      | 0.023                  | 0.050                 |
| 20                | 16.36                      | 0.07                                | 0.043                      | 0.028                  | 0.062                 |
| 22                | 19.80                      | 0.08                                | 0.051                      | 0.034                  | 0.075                 |
| 24                | 23.56                      | 0.09                                | 0.061                      | 0.040                  | 0.090                 |
| 26                | 27.65                      | 0.11                                | 0.072                      | 0.047                  | 0.105                 |
| 28                | 32.07                      | 0.13                                | 0.083                      | 0.055                  | 0.122                 |
| 30                | 36.82                      | 0.15                                | 0.096                      | 0.063                  | 0.140                 |
| 32                | 41.89                      | 0.17                                | 0.109                      | 0.071                  | 0.159                 |
| 34                | 47.29                      | 0.19                                | 0.123                      | 0.080                  | 0.180                 |
| 36                | 53.01                      | 0.21                                | 0.138                      | 0.090                  | 0.201                 |
| 40                | 65.45                      | 0.26                                | 0.170                      | 0.111                  | 0.249                 |
| 46                | 86.56                      | 0.35                                | 0.225                      | 0.147                  | 0.329                 |

Note: Amounts based on application of 200-ppm chlorine concentration into well water with Alkalinity of 100-ppm.

If well-water alkalinity, or the desired chlorine dosage level are greater than the standard values in the above table, adjust the amount of Chlorine Enhancer and the amount of hypochlorite concentrations as shown below.

| NuWell 410 Chlorine Enhancer | Hypochlorite                               |
|------------------------------|--|
| Amount above x (Alk/100)     | Amount above x (desired concentration/200) |

- Step 1: Determine Static Height (TD - SWL)
  - Step 2: Multiply Amt Chlorine Enhancer from table x Static Ht
  - Step 3: Determine Amount of hypochlorite product from table above
  - Step 4: Surface Batch = Water + NuWell 410 Chlorine Enhancer + Hypochlorite
- Example: Disinfect a 16 inch well, TD = 300 ft, SWL = 50 ft with Calcium Hypochlorite 65% active

- Step 1: Static Ht = (300ft - 50ft) = 250 ft
- Step 2: Amt Enhancer = 250 ft x 0.042 qts/ft = 11 quarts
- Step 3: Amt Hyp = 250 ft x 0.027 lbs/ft = 6.8 lbs.
- Step 4: Batch: 250 ft x 10.47 gal/ft = 2,618 gal water (2,618 gal + 11 quarts + 6.8 lbs)

TIP: Optimum results are obtained when the surface solution is 2-4 times the well volume, providing sufficient hypochlorite-ions to disperse into the gravel pack and immediate surrounding formation where coliform organisms exist. (For large wells this may not be practical and multiple batches with proportions of chemistry may be required to achieve the desired results.)

